


# June Practice Meeting Pack

- Week 1: Compression Hosiery Formulary ..... 1
- Week 2: PEFR Diary Template Improvement ..... 2
- Week 3: LTC Template ..... 3
- Week 4: Equipment Templates..... 4

## Week 1: Compression Hosiery Formulary

We've recently updated the compression hosiery and dressing formulary. Both of these can be accessed from the main 'Wound Management' template.



**Compression Hosiery Formulary**

Home Activa Actilymph Juxta Jobst Mediven Venotrain

ardens help & feedback

**Assessment:** Pulses + ABPI / Doppler

- $ABPI < 0.5$  = Compression stockings not to be worn
- $0.5-0.8$  = Only apply Class 1 compression
- $0.8-1.3$  = Compression stockings safe to wear
- $> 1.3$  = Avoid compression stockings

**Options:** Activa, Actilymph, Juxta, Jobst, Mediven, Venotrain

**Measurements:**

Left leg length: [ ] cm

mid thigh circumference: [ ] cm

calf circumference: [ ] cm

ankle circumference: [ ] cm

foot length: [ ]

Right leg length: [ ] cm

mid thigh circumference: [ ] cm

calf circumference: [ ] cm

ankle circumference: [ ] cm

foot length: [ ]

CKS

Current Hosiery: \_\_\_\_\_

Previous Hosiery: \_\_\_\_\_

New Acute, Drug Review, Allergies, Print Prescri..., Quick Print ...

Show recordings from other templates

Show empty recordings

Information, Print, Suspend, Ok, Cancel, Show Incomplete Fields

# June Practice Meeting Pack

## Week 2. PEFR Diary Template Improvement

The PEFR Diary template has been improved and when you print it out it will appear all on one page. You can access this from the Respiratory tab of the General Consultation template and also from the Asthma and COPD templates.

PEFR Diary
PEFR Diary used

		Peak Flow Rate												
		0	50	100	150	200	250	300	350	400	450	500	550	600
Day 1	am pm													
Day 2	am pm													
Day 3	am pm													
Day 4	am pm													
Day 5	am pm													
Day 6	am pm													
Day 7	am pm													
Day 8	am pm													
Day 9	am pm													
Day 10	am pm													
Day	am													

**Peak Expiratory Flow Rate Diary**

Do 3 Peak flow measurements in the morning and record the highest. Repeat in the

Information    Print    Suspend    Ok    Cancel    Show Incomplete Fields

PEFR Diary used

Date ▾    Checked    ...

No previous values

Show recordings from other templates


Show empty recordings

# June Practice Meeting Pack

## Week 3. LTC Template

We have updated the LTC Template. There is now a new tick box on the booking tab called 'Declined Invitation + Review'. Ticking this will stop the patient being sent 2nd and 3rd Invites.

We are also in discussions with WebPost and DocMail with regards to sending invitations out in bulk.



The screenshot shows the 'Long Term Condition Review' window with the 'Booking' tab selected. The main content area is titled 'Booking Appointment' and includes the following sections:

- Next Annual Review**: Due in June (31)
- Tests Required (+ not done in last 4m)**:
  - Blood tests required
  - COPD questionnaire required
  - ECG required
  - Spirometry required
  - Urine tests required
- Review Appointment**: General Review Required - Please book 30 minute appointment
- Conditions**:
  - Has Cardiovascular Conditions
  - Has Diabetes
  - Has Respiratory Conditions
- Pending Appointments**: (Empty list)

On the right-hand side, there is a sidebar with the following options:

- Free Slot Presets Search (31)
- Free Slot Search (31)
- Appointment Ledger
- New LTC Review task

A red button labeled 'Declined invitation + review' is visible in the sidebar, with a sub-note: '(stops patient being sent further invites)'. Below this, there are two checkboxes: 'Show recordings from other templates' (checked) and 'Show empty recordings' (unchecked).


At the bottom of the window, there are buttons for 'Information', 'Print', 'Suspend', 'Ok', 'Cancel', and 'Show Incomplete Fields'.

# June Practice Meeting Pack

## Week 4. Equipment Templates

The equipment templates have been updated and on them is a link to the 'Equipment Overdue' template. This allows you to easily send the patient an email or SMS reminder if their equipment is overdue.

This can be accessed via the 'Reception & Admin' template or via the 'BP Monitoring ', 'ECG Monitor' or 'Nebuliser' templates.



**Equipment Return Overdue**

ardens  
help & feedback

**Blood Pressure Monitor Needs Returning**  
We recently loaned a blood pressure monitor to you and it's return is now overdue.  
Please therefore return it as soon as possible to the surgery.  
If you have any questions, please contact us directly.

**ECG Monitor Needs Returning**  
We recently loaned an ECG monitor to you and it's return is now overdue.  
Please therefore return it as soon as possible to the surgery.  
If you have any questions, please contact us directly.

**Nebuliser Needs Returning**  
We recently loaned a nebuliser to you and it's return is now overdue.  
Please therefore return it as soon as possible to the surgery.  
If you have any questions, please contact us directly.

Launch Email  
 Launch SMS  
 Write Letter  
 Patient Contact Det...

To send by Email or SMS, highlight above > Copy with Ctrl+C  
Then click Launch Email/SMS > Paste with Ctrl + V > Send

Show recordings from other templates  
 Show empty recordings

Information Print Suspend Ok Cancel Show Incomplete Fields